



Recombinant Haemophilus influenzae Probable D-methionine transport system permease protein metI (metI)

Product Code	CSB-CF342005HTA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P46492
Product Type	Transmembrane Protein
Immunogen Species	Haemophilus influenzae (strain ATCC 51907 / DSM 11121 / KW20 / Rd)
Sequence	MWGVVATATYETVYISFASTLLAVLVGVPVGIWTFLLTGKNEILQNNRTHFVLNTI INIGR SIPFIILLILLPVTFRFIVGTVLGTTAAIPLSICAMPFVARLTANALMEIPNGLTEAA Q AMGATKWQIVRKFYLSEALPTLINGVTLLVTLVGVYSAMAGTQGGGGLGSLAIN YGRIRN MPYVTWVATIIIVLFVMISQKLGDTLAKKVDHR
Source	in vitro E.coli expression system
Target Names	metI
Protein Names	Recommended name: Probable D-methionine transport system permease protein metI
Expression Region	1-213
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.