



Recombinant Arabidopsis thaliana PRA1 family protein B2 (PRA1B2)

Product Code	CSB-CF879756DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9SIY7
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MSSSPAILPVTNQQAATQSQPPINSHAFRTFLSRLSSSLRESLSQRRPWLELV DRSSFAR PDSLTDSEFSRIRKNLAYFKVNYSIVSLVLAFLSLLSHPFSLVLLSLLGSWMFLY LFRSS DQPLVLFGRSFSRETLGLVLTIVVFMSTVGSLLTSALTIGIAIVCLHGAFRV PDDL FLDEQEPANAGLLSFIGNSAATSAAASVVAGRV
Source	in vitro E.coli expression system
Target Names	PRA1B2
Protein Names	Recommended name: PRA1 family protein B2 Short name= AtPRA1.B2
Expression Region	1-213
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.