



Recombinant Human Vesicle transport protein SFT2C (SFT2D3)

Product Code	CSB-CF682329HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q587I9
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MADLHRQLQEYLAQGGKAGGPAAAEPLLA AEKAE EPGDRPAEEWLGRAGLRW TWARSPAES AAAGLTCLPSVTRGQRLAAGGGCLLLAALCFGLAALYAPVLLLRARKFALLWSL GSALAL AGSALLRGG AACGRLLRCEEAPSRPALLYMAALGATLFAALGLRSTLLTVLGA GAQVAAL LAALVGLLPWGGGTALRLALGRLGRGAGLAKVLPV
Source	in vitro E.coli expression system
Target Names	SFT2D3
Protein Names	Recommended name: Vesicle transport protein SFT2C Alternative name(s): SFT2 domain-containing protein 3
Expression Region	1-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.