



Recombinant Schizosaccharomyces pombe Secretory component protein psh3 (psh3)

Product Code	CSB-CF896973SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9Y876
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MAKRSIFRFADEKGLKVAARYGVLMSTSFIFALLFHSSVADVNTLWSPGPESAF DAAETY YTLVAGSHFIVKYTVYTIMGLNMIFHLIQATGAKGDDKLVFFYSSTLLYLTLALILFIV NVA PSMLVVKLQNYVQFPRNMHLSVLAASHVLVEFLLAGVILIQLGYVFGYHVQSIQ QREYAE DMREQELAEKAKLESESATTQSVETVSTESVSKRK
Source	in vitro E.coli expression system
Target Names	psh3
Protein Names	Recommended name: Secretory component protein psh3
Expression Region	1-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.