



Recombinant Danio rerio Claudin-7-B (cldn7b)

Product Code	CSB-CF897044DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9YH92
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MAHKGLQLLGFTLSLLGLIGLIIGTIMPQWKMSAYVGDNIITAIAMYQGLWMSCA YQSTG QQQCKVYDSVLQLDSALQATRALMVVAILLTVAGLGVASMGMKCTNCGGDDK VKKSRIAM TGGIILSVGALCSIVACGWFTSQIIRDFYNPFTPVNTKYEFGAAIFIAWAGAFLDI MGGG MLASSCSKQSSPNYPKSSRPVKSSRPPSSSKEYV
Source	in vitro E.coli expression system
Target Names	cldn7b
Protein Names	Recommended name: Claudin-7-B Alternative name(s): Claudin-like protein ZF4A22
Expression Region	1-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.