



Recombinant Arabidopsis thaliana PRA1 family protein B6 (PRA1B6)

Product Code	CSB-CF868075DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9LYQ4
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MASPLLPTSTTPDQLPGGDPQLLSSLRVLLSRVLATVRHASADARPWAELVDR SAFSRPP SLSEATSRVRKNFSYFRANYITLVAILLAASLLTHPFALFLLASLAASWLFYFFR PADQ PLVIGGRTFSDLETLGILCLSTVVVMFMTSVGSLLMSTLAVGIMGVAIHGAFRAP EDLFL EEQEAIGSGLFAFFNNNASNAAAAAIATSMSRVRV
Source	in vitro E.coli expression system
Target Names	PRA1B6
Protein Names	Recommended name: PRA1 family protein B6 Short name= AtPRA1.B6 Alternative name(s): Prenylated Rab acceptor 3
Expression Region	1-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.