



Recombinant *Schizosaccharomyces pombe* Uncharacterized membrane protein C30D10.09c (SPBC30D10.09c)

Product Code	CSB-CF522679SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O14355
Product Type	Transmembrane Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MQLLMLEQISLASSVVATTVLVAPVLSTIINLLTNFGQRIFIAADSSKSTSMEFLS TIIG AGYPIYKTYLLELPSKRSQLPKAFQLRNEEHKSIEEERRRLMAYWCVYGCVT AAESIL GRFLSWVPFYSTSKIVFWLWLLNPRTQGAAFIYASYISPFLSDHKAANNFLEKL VQFTT RQPLVLNAWALVKSLIDKLPKGDVEAPGSDADTKKSK
Source	in vitro E.coli expression system
Target Names	SPBC30D10.09c
Protein Names	Recommended name: Uncharacterized membrane protein C30D10.09c
Expression Region	1-217
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.