



Recombinant Arabidopsis thaliana PRA1 family protein B3 (PRA1B3)

Product Code	CSB-CF880783DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9FLB6
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MMANPPTLPISDHSGGGSQSQQPVSTPAFRTFLSRLSSSIRQSLSQRRPWLEL VDRSAIS RPESLTDAYSRI RRNLPHYFKVNYVTIVSLVLALSLLSHPFLLVLLCLFCAWIFLY LFRP SDQPLVVLGRTFSDRETLGVLVILTIVVVFLTSVGSLLTSALMIGFGIVCLHGAFR VPED LFLDDQEPANTGLLSFLSGAATSAAVAAASTPASGRV
Source	in vitro E.coli expression system
Target Names	PRA1B3
Protein Names	Recommended name: PRA1 family protein B3 Short name= AtPRA1.B3 Alternative name(s): Prenylated Rab acceptor 2
Expression Region	1-217
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.