



Recombinant Mouse Leukocyte surface antigen CD53 (Cd53)

Product Code	CSB-CF736818MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q61451
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	GMSSLKLLKYVLFIFNLLFWVCGCCILGFGIYFLVQNTYGVLFERNLPFLTGLNILV IVGS IIMVVAFLGCMGSIKENKCLLMSFFVLLLIILLAEVTIAILLFVYEQKLNTLVAEGLN DS IQHYHSDNSTMKAWDFIQTQLQCCGVNGSSDWTSGPPSSCPSGADVQGCYN KAKSWFHSN FLYIGIITICVCVIQVLGMSFALTLCNCQIDKTSQALGL
Source	in vitro E.coli expression system
Target Names	Cd53
Protein Names	Recommended name: Leukocyte surface antigen CD53 Alternative name(s): Cell surface glycoprotein CD53 CD_antigen= CD53
Expression Region	2-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.