



Recombinant Inner membrane protein YghB (yghB)

Product Code	CSB-CF364441EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AA62
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MAVIQDIIAALWQHDFaalADPHIVSVVYFVMFATLFLenGLLPASFLPGDSSLIL AGAL IAQGVMDFLPTIAILTAASLGCWLSYIQGRWLGNTKTKVKGWLAQLPAKYHQR ATCMFDR HGLLALLAGRFLAFVRTLLPTMAGISGLPNRRFQFFNWLSGLLWVSVVTSFGY ALSMIPF VKRHEDQVMTFLMILPIALLTAGLLGTLFVVIKKKYCNA
Source	in vitro E.coli expression system
Target Names	yghB
Protein Names	Recommended name: Inner membrane protein YghB
Expression Region	1-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.