



Recombinant Danio rerio Transmembrane protein 179B (tmem179b)

Product Code	CSB-CF023762DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7T392
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MALPRPLLELLLYGSCFICGIITAASVTISQGSFDGKCMLYGSVRLNSSTIDVLS FSSL SLCYFVSAISVCVAVFCFSLTLYWIYIAFVDGEIKREKLWMNVTLGLSGVFLFFL LVTGC ILKIGRDRLCDSLLHTVSNITRCEEAAQNKSWRSPINASQFYTRLHSAETA VVWN FFFWMI IVVLVLIQRHKGSEIRPGTEDPSAPPSETEPFFNRPGRP
Source	in vitro E.coli expression system
Target Names	tmem179b
Protein Names	Recommended name: Transmembrane protein 179B
Expression Region	1-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.