



Recombinant Danio rerio Transmembrane protein 17A (tmem17a)

Product Code	CSB-CF703475DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q502E0
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MPVFYTPIPQNLRVGLANVSGSVFINNRTDRDSADFKDHYLENEACAEVSSHLP LQMMLYF NMFFFPFWWISSELLMLQLKFSYLPVYYQCLLVGTGMVLISIFEVLRMYLGYAGNL KEKVPE LAGFWLISFLFQLPILLFFITDPDIILPLERAVHSLYLAFLLGELMASFLALRVMTR KL AQQFHMRQFGHVQGLHTAEALPMFGLPYGGRSVLPVQDI
Source	in vitro E.coli expression system
Target Names	tmem17a
Protein Names	Recommended name: Transmembrane protein 17A
Expression Region	1-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.