



Recombinant Rat Claudin-3 (Cldn3)

Product Code	CSB-CF723578RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q63400
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MSMGLEITGTSLAVLGWLVCTIVCCALPMWRVSAFIGSSIITAQITWEGLWMNCV VQSTGQ MQCKMYDSSLALPQDLQAARALIVVSILLAAGLLVALVGAQCTNCVQDETA KITIVA GVLFLAAVLTLPVSVSANTIIIRDFYNPLVPEAQKREMGTGLYVGWAAAALQ LLGGALL CCSCPPREKYAPTILYSAPRSTGPGTGTGTAYDRKDYY
Source	in vitro E.coli expression system
Target Names	Cldn3
Protein Names	Recommended name: Claudin-3 Alternative name(s): Rat ventral prostate.1 protein Short name= RVP1
Expression Region	1-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.