



Recombinant Human Leukocyte surface antigen CD53 (CD53)

Product Code	CSB-CF004944HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P19397
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MGMSSLKLLKYVLFNLLFWICGCCILGFGIYLLIHNNFGVLFHNLPSLT LGNV FVIVG SIIMVVAFLGCMGSIKENKCLLMSFFILLIILLAEVTLAILLFVYEQKLNEYVAKGL TD SIHRYHSDNSTKAAWDSIQSFLQCCGINGTSDWTSGPPASCPSDRKVEGCYA KARLWFHS NFLYIGIITICVCVIEVLGMSFALT LNCQIDKTSQTIGL
Source	in vitro E.coli expression system
Target Names	CD53
Protein Names	Recommended name: Leukocyte surface antigen CD53 Alternative name(s): Cell surface glycoprotein CD53 Tetraspanin-25 Short name= Tspan-25 CD_antigen= CD53
Expression Region	1-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.