



Recombinant Inner membrane protein YqjA (yqjA)

Product Code	CSB-CF359407EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AA64
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Sequence	MELLTQLLQALWAQDFETLANPSMIGMLYFVLFVILFLENGLLPAAFLPGDSSLV LVGVL IAKGAMGYPQTILLTVAASLGCWVSYIQGRWLGNTTRTVQNWLSHLPAHYHQ RAHHLFHK HGSLALLIGRFIAFVRTLLPTIAGLSGLNNAFQFFNWMSGLLWVILTTLGMYL GKTPV FLKYEDQLMSCLMLLPVLLVFLAGSLVVLWKKKYGNRG
Source	in vitro E.coli expression system
Target Names	yqjA
Protein Names	Recommended name: Inner membrane protein YqjA
Expression Region	1-220
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.