



Recombinant *Saccharomyces cerevisiae* Protein ORM1 (ORM1)

Product Code	CSB-CF345664SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53224
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MTELDYQGTAEEAASSTYSRNQTDLKPFPSPAGSASSSIKTTEPVKDHRRRRSSS IISHVEP ETFEDENDQQLLPNMNATWVDQRGAWIIHVVIILLKLFYNLFPGVTTTEWSWTL TNMTYV IGSYVMFHLLIKGTPFDNFNGGAYDNLTMWEQIDDETLYTPSRKFLISVPIALFLVS THYAH YDLKLFSWNCFLTTFGAVVPKLPVTHRLRISIPGITGRAQIS
Source	in vitro E.coli expression system
Target Names	ORM1
Protein Names	Recommended name: Protein ORM1
Expression Region	1-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.