



Recombinant Danio rerio Calcium release-activated calcium channel protein 1 (orai1)

Product Code	CSB-CF761459DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6TLE6
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MKMSRSEHSLQALSWRKLYLSRAKLKASSRTSALLSGFAMVAMVEVQLDTNH DYPPGLLI AFSACTTVLVAVHLFALMVSTCILPNIEAVSNVHNLNSVKESPHERMHRHIELA WAFSTV IGTLLFLAGVVLLCWVKFLPIKPKGQKNGTVSAGVAAAITSTSIMVPFGLVFIVFA VHFY RSLVSHKTDRQFQELEELEDLQNELDHREEVSTLQSPGSLYP
Source	in vitro E.coli expression system
Target Names	orai1
Protein Names	Recommended name: Calcium release-activated calcium channel protein 1 Alternative name(s): Protein orai-1 Transmembrane protein 142A
Expression Region	1-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.