



Recombinant Putative protein-disulfide oxidoreductase (Ecok1_30590)

Product Code	CSB-CF372343EJE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A1AFW3
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O1:K1 / APEC
Sequence	MGIKGMWKDLRTSPVDLVRWQEQRLLWLLMAVAMGALIILAHSFFQIYLYMA PCEQCVY IRYAMFVMVIGGLVAAINPKNIILKLIGCVMAFYGSILGLKFSCLKNDIHHAVHNP DPDS LFGVQGCSTDPTFPFNLPLAQWAPNWFKPTGDCGYDAPIVPDGVTLSSSTQQ WFVEMYQQS EGWYLLPPWHFMNMAQACMLAFGMCLVLLVIMSGAWALKIIRG
Source	in vitro E.coli expression system
Target Names	dsbl
Protein Names	Recommended name: Putative protein-disulfide oxidoreductase
Expression Region	1-223
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.