



Recombinant Human Gap junction beta-7 protein (GJB7)

Product Code	CSB-CF757637HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6PEY0
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MSWMFLRDLLSGVKNKYSTGTGWIWLAVVFVFRLLVYMVAAEHVWKDEQKEF ECNSRQPGC KNVCFDDFFPISQVRLWALQLIMVSTPSLLVLHVAYHEGREKRHRKKLYVSP GTMDGGL WYAYLISLIVKTGFEIGFLVLFYKLYDGFSVPYLIKCDLKPCPNTVDCFISKPTK TIFI LFLVITSCLCIVLNFIELSFLVLKCFIKCCLQKYLKKPQVLSV
Source	in vitro E.coli expression system
Target Names	GJB7
Protein Names	Recommended name: Gap junction beta-7 protein Alternative name(s): Connexin-25 Short name= Cx25
Expression Region	1-223
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.