



Recombinant *Drosophila melanogaster* Uncharacterized FAM18-like protein CG5021 (CG5021)

Product Code	CSB-CF815534DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8IQC1
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	MASATVRNVPLDDDDTIPFGEEDEMRDPSRAGQKYTHPYVTFFHLFFRGAAILI YMF CGW FSDSFITSFV FVVLFLSADFWTVKNISGRLLVGLRWVNYVDDDGVSHWVFESK NSESYSQS RVNKNEQRIFWLGLILCPVFWGLFFLFALFGLKFKWLLLVMIAIALNAANLYGYV KCNYG ANKDLNSAATDFVKTQFFKNAVDIMTRPSGAPPPTNVRPTGVV
Source	in vitro <i>E.coli</i> expression system
Target Names	CG5021
Protein Names	Recommended name: Uncharacterized FAM18-like protein CG5021
Expression Region	1-223
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.