



Recombinant Human Synaptogyrin-2 (SYNGR2)

Product Code	CSB-CF023012HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O43760
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MESGAYGAAKAGGSFDLRRFLTQPQVVARAVCLVFALIVFSCIYGEGYSNAHE SKQMYCV FNRNEDACRYGSAIGVLAFLASAFFLVVDAYFPQISNATDRKYLVIDLLFSAL WTFLWF VGFCFLTQWAVTNP KDVLVGADSVRAAITFSFFSIFSWGVLASLAYQRYKAG VDDFIQN YVDPTPDPNTAYASYPGASVDNYQQPPFTQNAETTEGYQPPPVY
Source	in vitro E.coli expression system
Target Names	SYNGR2
Protein Names	Recommended name: Synaptogyrin-2 Alternative name(s): Cellugyrin
Expression Region	1-224
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.