



Recombinant Cat CD9 antigen (CD9)

Product Code	CSB-CF004969CA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P40239
Product Type	Transmembrane Protein
Immunogen Species	Felis catus (Cat) (Felis silvestris catus)
Sequence	PVKGGTKCIKYLLFGFNFIWLAGIAVLAVGLWLRFDSQTKSIFEQDSQPSSFY TGVYIL IGAGALMMLVGFLGCCGAVQESQCMLGLFFGFLLVIFAIEIAAAIWGYSHKDEVI QEVQE FYKDTYNKLKSKDEPQRDTLKAHYALDCCGLAGGVEQFISDICPQKDILSSITV KPCPE AIKEVFHNKFHIIIGAVGIGIAVVMIFGMIFSMILCCAIRRSREMV
Source	in vitro E.coli expression system
Target Names	CD9
Protein Names	Recommended name: CD9 antigen Alternative name(s): CD_antigen= CD9
Expression Region	2-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.