



Recombinant Escherichia coli Respiratory nitrate reductase 1 gamma chain (narI)

Product Code	CSB-CF318854ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P11350
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	<pre>MQFLNMFFFDIYPYIAGAVFLIGSWLRYDYGQYTWRAASSQMLDRKGMNLAS NLFHIGIL GIFVGHHFFGMLTPHWMYEAWLPIEVKQKMAMFAGGASGVLCCLIGGVLLLKRRL FSPRVRA TTTGADILILSLLVIQCALGLLTIPFSAQHMDGSEMMKLVGWAQSVVTFHGGAS QHLDGV AFIFRLHLVLGMTLFLFPFSRLIHIWVSPVEYLTRKYQLVRARH</pre>
Source	in vitro E.coli expression system
Target Names	narI
Protein Names	Recommended name: Respiratory nitrate reductase 1 gamma chain EC=1.7.99.4 Alternative name(s): Cytochrome B-NR Nitrate reductase A subunit gamma Quinol-nitrate oxidoreductase subunit gamma
Expression Region	1-225
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.