



Recombinant Mouse Clarin-3 (Cln3)

Product Code	CSB-CF805895MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8BHH8
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MPTTQKTLMFLSGFLTSLGSVVICILATQAWITSRIFFTDAISNGTIVITYGLFR GTS AQELNEGLQDLKDNFEVLGILDNSSQKSLHLVVILLILSLAASVLSSVFTFYNSI SNPY QTFLGPMGVYTWNGLSASFVFLAMVLFVGNNAESNHLSDKLSQKLYPDTTNKR TTHTYGYF FWLTLHVIFLNIVTAVIIIFYQKARYRQKQEQRKPVEYAPRDGILF
Source	in vitro E.coli expression system
Target Names	Cln3
Protein Names	Recommended name: Clarin-3 Alternative name(s): Transmembrane protein 12
Expression Region	1-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.