



Recombinant Thioalkalivibrio sp. Electron transport complex protein RnfE (rnfE)

Product Code	CSB-CF490631TNX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B8GMB9
Product Type	Transmembrane Protein
Immunogen Species	Thioalkalivibrio sulfidophilus (strain HL-EbGR7)
Sequence	MSDISYREINANGFWHNNPGLVQLLGLCPLLAIISGTVVNALGLGLATTLTLVAS NVTVSL IRHWVRPEIRIPVFLIIASVVTAIELAMNAFFHELYLILGIFIPLIVTNCAIIGRAEAF ASKQPIPKALADGLAMGLGFTCVLVALGALREAVGHGTLLADAHLMFGEAARG FSLTLFE EYRGFLLALLPPGAFIALGLLLIALKNIIDARLQKRQPAQAAVPVEAAG
Source	in vitro E.coli expression system
Target Names	Tgr7_2631
Protein Names	Recommended name: Electron transport complex protein RnfE
Expression Region	1-228
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.