



Recombinant Acorus calamus Chloroplast envelope membrane protein (cemA)

Product Code	CSB-CF661092ACE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q3V523
Product Type	Transmembrane Protein
Immunogen Species	Acorus calamus (Sweet flag)
Sequence	MPKKKGFTPLPYLASIVFLPWWVSLSFNKSLEPWVTNWWNTRQSETFLNDIQ ERNVLERF IELEELFLLDEMLKENPETRMKNLRIGIHNETIQLVKTDNEYHLHTILHFSTNIICF AIL SVYSILGNEELVILNSWVQEFLYNLSDTIKAFSILLVTDLWIGFHSPHGWELMIG SVYND FGLAHNEQIISGLVSTFPVILDITIVKYWIFHYLNRVSPSLVVIYHSMND
Source	in vitro E.coli expression system
Target Names	cemA
Protein Names	Recommended name: Chloroplast envelope membrane protein
Expression Region	1-229
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.