



Recombinant Bacillus subtilis Uncharacterized protein yknW (yknW)

Product Code	CSB-CF516298BRJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O31709
Product Type	Transmembrane Protein
Immunogen Species	Bacillus subtilis (strain 168)
Sequence	METNVEKNSGTATEKPSLFGVITSPSVQFERIRERPAVWGPLLIVAAIIIVGAVL QSLGT DYSELLKSQDTQGLSAEQMETVATITKFGGMAGAIIGGIAALFIAPLIYWLCVKV SGGVT TYKKMLSLSLFFVSLISSLGLLVNGIVAFTTVDVNPLYSTTSLAGIIPSDGALASVLN TFEI FSIWSFVLLAIGLHKTGGISKKAGWISAILFGILVVFLFSGLINSVAGA
Source	in vitro E.coli expression system
Target Names	yknW
Protein Names	Recommended name: Uncharacterized protein yknW
Expression Region	1-231
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.