



Recombinant Danio rerio Transmembrane protein 179 (tmem179)

Product Code	CSB-CF023761DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6AZD1
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MAVDNFLFGQCILYFLAFLFGFIAVVPLSENGDDFQGKCLLFTEGIWQENENMT MGKQRFI VEEWGPESSECRFITFVGIVSLILSAVQAWRTFFFFLCKGHDDSLFHSFLNLLLSLL VLFVV FVAGTISSVGFSIWCDSVTENGAMPSSCEDLQDLDLELGVENSSFYDQFAIAQ FGLWSAW LCWLGLTVLAFVKVYHNHRQQELLESLVQEKELLLGHPLQRSSYVYNRNAMI
Source	in vitro E.coli expression system
Target Names	tmem179
Protein Names	Recommended name: Transmembrane protein 179
Expression Region	1-232
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.