



Recombinant Mouse Lysosomal-associated transmembrane protein 4A (Laptm4a)

Product Code	CSB-CF726690MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q60961
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MVSMTFKRSRSDRFYSTRCCGCFHVRTGTIILGTWYMVVNLLMAILLTVETH PNSMPAV NIQYEVIGNYYSSERMADNACVLFVSVLMFISSMLVYGAIQVGVWLPFFCY RLDFD VLSCLVAISLTYLPRIKEYLDQLPDFPYKDDLALDSSCLLFIVLVFFVVFIIFKAY LI NCVWNCYKYINNRNVPEIAVYPAFETPPQYVLPTYEMAVKIPEKEPPPPYLP
Source	in vitro E.coli expression system
Target Names	Laptm4a
Protein Names	Recommended name: Lysosomal-associated transmembrane protein 4A Alternative name(s): Golgi 4-transmembrane-spanning transporter Mouse transporter protein Short name= MTP
Expression Region	1-233
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.