



Recombinant *Oryza sativa* subsp. *indica* CASP-like protein Osl_27482 (Osl_27482)

Product Code	CSB-CF383743OFF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2YQB8
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Sequence	<p>MSGSDTSGSVHVDEHGHGKASSSYDGAGAPAPAPAPFQGHRKAGSGSSDV PFLLRSGGSG GDGLRRCLGLIDFVLRVAAFPGPTLAAAISIGTSDERLSVFTNYFQFRARFDDFP AFEFFI VANAIAAGYMLSLPFS AATIMSSKATGVKLLLLICDTIMVGLLTAAASAAAAMV YVAHE GNLRANWVPICLQFHGFCQRTSGAVIASFLAVFVLMVLIVMAAFTMPRRTHHT AS</p>
Source	in vitro <i>E.coli</i> expression system
Target Names	Osl_27482
Protein Names	Recommended name: CASP-like protein Osl_27482
Expression Region	1-235
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>