



Recombinant Human Derlin-3 (DERL3)

Product Code	CSB-CF836283HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q96Q80
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAWQGLAAEFLQVPAVTRAYTAACVLTTAAVQLELLSPFQLYFNPHLVFRKFQ VWRLVTN FLFFGPLGFSSFFNMLFVFRYCRMLEEGSFRGRTADFVFMFLFGGVLMTLLGL LGSLFFL GQALMAMLVYVWSRRSPRVRVNFFGLLTFQAPFLPWALMGFSLLLGNSILVDL LGIAVGH IYYFLEDVFPNQPGGKRLQLTPGFLKLLLDAPAEDPNYLPLPEEQPGPHLPPPQ Q
Source	in vitro E.coli expression system
Target Names	DERL3
Protein Names	Recommended name: Derlin-3 Alternative name(s): Degradation in endoplasmic reticulum protein 3 Short name= DERtrin-3 Der1-like protein 3
Expression Region	1-235
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.