



Recombinant Rat CD81 antigen (Cd81)

Product Code	CSB-CF726808RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q62745
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MGVEGCTKCIKYLFFVFNFWLGGVILGVALWLRHDPQTTTTLLYLELGDKPA PSTFYV GIYILIAVGAVMMFVGFVGLGCGYGAIQESQCLLGTFFFTCLVILFACEVAAGIWGFVN KDQIA KDVKQFYDQALQQAVMDDDANNAKAVVKTFFHETLNCCGSNTLTTTAVLRN SLCPSSSN SFTQLLKEDCHQKIDELFSGKLYLIGIAAIVVAVIMIFEMILSMVLCCGIRNSSVY
Source	in vitro E.coli expression system
Target Names	Cd81
Protein Names	Recommended name: CD81 antigen Alternative name(s): 26 kDa cell surface protein TAPA-1 Target of the antiproliferative antibody 1 CD_antigen= CD81
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.