



Recombinant Pan troglodytes CD81 antigen (CD81)

Product Code	CSB-CF004960EQV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P60034
Product Type	Transmembrane Protein
Immunogen Species	Pan troglodytes (Chimpanzee)
Sequence	MGVEGCTKCIKYLLFVFNFWLWLAGGVILGVALWLRHDPQTTNLLYLELGDKP APNTFYV GIYILIAVGAVMMFVGFLGCGYGAIQESQCLLGTTFFTCVLVIFACEVAAGIWGFVN KDQIA KDVKQFYDQALQQAVVDDDANNAKAVVKTFFHETLDCCGSSTLTALTTSVLKN NLCPSGSN IISNLFKEDCHQKIDDLFSGKLYLIGIAAIVVAVIMIFEMILSMVLCCGIRNSSVY
Source	in vitro E.coli expression system
Target Names	CD81
Protein Names	Recommended name: CD81 antigen Alternative name(s): 26 kDa cell surface protein TAPA-1 Target of the antiproliferative antibody 1 CD_antigen= CD81
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.