



Recombinant Escherichia coli TVP38/TMEM64 family membrane protein ydjX (ydjX)

Product Code	CSB-CF304338ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76219
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MNAERKFLFACLIFALVIYAIHAFGLFDLLTDLPHLQTLIRQSGFFGYSLYILLFIIA TL LLPGSILVIAGGIVFGPLLGTLLSLIAATLASSCSFLLARWLGRDLLLLKYVGHSN TFQA IEKGIARNGIDFLILTRLIPLFPYNIQNYAYGLTTIAFWPYTLISALTTLPGIVIYTVM A SDLANEGITLRFILQLCLAGLALFILVQLAKLYARHKHVDLSASRRSPLTHPKNE G
Source	in vitro E.coli expression system
Target Names	ydjX
Protein Names	Recommended name: TVP38/TMEM64 family membrane protein ydjX
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.