



Recombinant *Oryza sativa* subsp. *japonica* CASP-like protein Os08g0101900 (Os08g0101900, LOC_Os08g01160)

Product Code	CSB-CF747751OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6Z1Y7
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Sequence	MSGSDTSGSVHVDHGHGHGKASSSYDGAGAPAPAPAPFQGHRKAGSGSS DVPFLLRSGG SGGDGLRRCLGLIDFVLRVAAFGPTLAAASIGTSDERLSVFTNYFQFRARFDD FPAFEF FIVANAIAAGYMLSLPFS AATIMSSKATGVKLLLLICDTIMVGLLTAAASAAAAM VYVA HEGNLRANWVPICLQFHGFCQRTSGAVIASFLAVFVLMVLIVMAAFTMPRRTH HTAS
Source	in vitro <i>E.coli</i> expression system
Target Names	Os08g0101900
Protein Names	Recommended name: CASP-like protein Os08g0101900
Expression Region	1-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.