



Recombinant Rat CD63 antigen (Cd63)

Product Code	CSB-CF004950RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P28648
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	AVEGGMKCVKFLLYVLLLAFCACAVGLIAIGVAVQVVLKQAITHETTAGSLLPVV IIAVG AFLFLVAFVVGCCGACKENYCLMITFAIFLSLIMLVEVAVAIAGYVFRDQVKSEFS KSFQK QMQNYLTDNKTATILDKLQKENKCCGASNYTDWERIPGMAKDRVPDSCCINIT VGCGNDF KESTIHTQGCVETIAAWLRKNVLLVAGAALGIAFVEVLGIIFSCCLVKSIRSGYEV M
Source	in vitro E.coli expression system
Target Names	Cd63
Protein Names	Recommended name: CD63 antigen Alternative name(s): Mast cell antigen AD1 CD_antigen= CD63
Expression Region	2-238
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.