



Recombinant Histidine transport system permease protein hisM (hisM)

Product Code	CSB-CF365200EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AEU5
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MIEILHEYWKPLLWTDGYRFTGVAILWLLILSVVIGGV LALFLAIGRVSSNKYIQ FPIW LFTYIFRGTPLYVQLLVFYSGMYTLEIVKGTFLNAFFRSGLNCTVLALTLNTCA YTTEI FAGAIRSVPHGEIEAARAYGFSTFKMYRCIILPSALRIALPAYSNEVILMLHSTAL AFTA TVPDLLKIARDINAATYQPFTAFGIAAVLYLIISYVLISLFRRRAEKRWLQHVKPSST H
Source	in vitro E.coli expression system
Target Names	hisM
Protein Names	Recommended name: Histidine transport system permease protein hisM
Expression Region	1-238
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.