



Recombinant Mouse Dolichyldiphosphatase 1 (Dolpp1)

Product Code	CSB-CF867562MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9JMF7
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MAADGQCCLPASWRPVTLTHVEYPAGDLSGHELLAYLSLSPIFVVVGFLTLIIFKR ELHTI SFLGGLALNQGVNWLKIHVIQEPKPCGGPHTAVGTKYGMPSHSHSQFMWFFSV YSFLFLYL RMHQTNNARFLDLLWRHVLSLGLLTAFLVSYSRVYLLYHTWSQVFGGVAG SLMAVAWF IITQEILTPLFPRIAAPWPISEFFLIRDTSIPNVLWFEYTVTRAEARNRQRKLGTKL Q
Source	in vitro E.coli expression system
Target Names	Dolpp1
Protein Names	Recommended name: Dolichyldiphosphatase 1 EC= 3.6.1.43 Alternative name(s): Dolichyl pyrophosphate phosphatase 1 Protein 2-23
Expression Region	1-238
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.