



Recombinant Dictyostelium discoideum Probable derlin-1 homolog (derl1)

Product Code	CSB-CF709711DKK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q54IC9
Product Type	Transmembrane Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Sequence	MDGVKEWFNSIPPVSRYMFAIFLGIPVLAAMHLISFNLYLDFTFTFKHFHLWRL ITAPC ISSLGPMFLFNLIFFYQYTTTRLESLNYAGKSDDYLCIIFISICNIIFGLIFEYYFLG T MTIMSLIYIYSRMNPTGTSNFYGFFSFKTIYLPWVFLVAHFLQTGHPPYSDFLAI VSGHI FFYLTDIYPRANGVPALIKTPKFITNIFNKGDRNPNNVRRDPRTGRPIQEGGYN WGQGHA LG
Source	in vitro E.coli expression system
Target Names	derl1
Protein Names	Recommended name: Probable derlin-1 homolog
Expression Region	1-242
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.