



Recombinant Human Sarcospan (SSPN)

Product Code	CSB-CF622786HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q14714
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MGKNKQPRGQQRQGGPPAADAAGPDDMEPKKGTGAPKECGEEEPRTCCG CRFLLLALLQ LALGIAVTVVGFLMASISSLLVRDTPFWAGIIVCLVAYLGLFMLCVSYQVDERT CIQFS MKLLYFLLSALGLTVCVLAVAFAAHHYSQLTQFTCETTLDSCQCKLPSSSEPLSR TFVYRD VTDCTSVTGTGFKLFLLIQMILNLVCGLVCLLACFVMWKHRYQVFYVGVRICSLT ASEGPQ QKI
Source	in vitro E.coli expression system
Target Names	SSPN
Protein Names	Recommended name: Sarcospan Alternative name(s): K-ras oncogene-associated protein Kirsten-ras-associated protein
Expression Region	1-243
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.