



Recombinant Zea mays Derlin-1.1 (DER1.1)

Product Code	CSB-CF671208ZAX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q4G2J6
Product Type	Transmembrane Protein
Immunogen Species	Zea mays (Maize)
Sequence	MSSPAEYYKSLPPIISKAYGTLCCFFTTVLVQLQILHPLFLYLDYPLVFKKFEIWRL TSFF FLAPFSMKFGIRLLMIARYGVMLEKGAFDKRTADFLWMMIFGAISLLVLSIPLFN SFFL GIPMVSMMLLYVWSRENPNQAQINIYGLVQLRSFYLPWAMLLLDVIFGSSLMPGLL GIMVGH LYYFFAVLHPLATGKSYLKTPKWVHKIVARFRIGMQANSPVRPPANGNSGSGV FRGRSYR LNQ
Source	in vitro E.coli expression system
Target Names	DER1.1
Protein Names	Recommended name: Derlin-1.1 Alternative name(s): ZmDerlin1-1
Expression Region	1-243
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.