



Recombinant *Gossypium barbadense* ATP synthase subunit a, chloroplastic (atpl)

Product Code	CSB-CF372057GEP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A0ZZ23
Product Type	Transmembrane Protein
Immunogen Species	<i>Gossypium barbadense</i> (Sea-island cotton) (Egyptian cotton)
Sequence	MNGISNALNGLYDISGVEVGQHFYWQIAGFQVHAQVLITSWVIAILLGSAVI AVRNPQT IPTAGQNFFFEYVLEFIRDVSKTQIGEEYGPWVVPFIGTMFLFIFVSNWSGALLPW KIIQLP HGELAAPTNDINTTVALALLTSVAYFYAGLSKKGLGYFSKYIQPTPILLPINILED FTKP LSLSFRLFGNILADELVVVVLVSLVPSVVPVPMFLGLFISGIQALIFATLAAAYIG ESM EGH
Source	in vitro E.coli expression system
Target Names	atpl
Protein Names	Recommended name: ATP synthase subunit a, chloroplastic Alternative name(s): ATP synthase F0 sector subunit a F-ATPase subunit IV
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.