



Recombinant Human Claudin-12 (CLDN12)

Product Code	CSB-CF005493HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P56749
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MGCRDVHAATVLSFLCGIASVAGLFAGTLLPNWRKLRRLITFNRNEKNLTVYTGL WVKCAR YDGSSDCLMYDTTWYSSVDQLDLRVLQFALPLSMLIAMGALLLCLIGMCNTAF RSSVPNI KLA KCLVNSAGCHLVAGLLFFLAGTVSLSPSIWVIFYNIHLNKKFEPVFSFDYAV YVTIA SAGGLFMTSLILFIWYCTCKSLPSPFWQPLYSHPPSMHTYSQPYSARSRLSAIE IDIPVV SHTT
Source	in vitro E.coli expression system
Target Names	CLDN12
Protein Names	Recommended name: Claudin-12
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.