



Recombinant Derlin-1 (cup-2)

Product Code	CSB-CF842369CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q93561
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MDLENFLLGIPIVTRYWFLASTIIPLLGRFGFINVQWMFLQWDLVVNKFQFWRP LTALIY YPVTPQTGFHWLMMCYFLYNYSKALESETYRGRSADYLFMLIFNWWFFCSGLC MALDIYFL LEPMVISVLYVWCQVNKDTIVSFWFGMRFPARYLPWVLWGFNAVLRGGGTNE LVGILVGH AYFFVALKYPDEYGVDLISTPEFLHRLIPDEDGGIHGQDGNIRGARQQPRGHQ WPGGVGA RLGGN
Source	in vitro E.coli expression system
Target Names	cup-2
Protein Names	Recommended name: Derlin-1 Alternative name(s): Coelomocyte uptake defective protein 2 DER1-like protein 1 cDerlin-1
Expression Region	1-245
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.