



Recombinant *Drosophila melanogaster* Derlin-1 (Der-1)

Product Code	CSB-CF894261DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VQ57
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	MDAGVWYRSLPRFTRYWLTATVVL SMLCRFDVIPLHWHLDRSAVFSKLQLW RCMTSLFV FPISSNTAFHFLINCFIVQYSSKLEKDQYSRSPADYLYLLIVSAVLANIGGMIFN VYFL MDTLVLAITYIWCQLNKDVTVSFWFGTRFKAMYLPWVLA AFEFIFHFSLASLVGI FVGHV YYFFKFQYSQDLGGTLPLETPQFLKRLVPDVSGGFGGFGGLPPESRAPPRQAT ESPWGRGM TLGRN
Source	in vitro <i>E.coli</i> expression system
Target Names	Der-1
Protein Names	Recommended name: Derlin-1 Alternative name(s): DER1-like protein 1
Expression Region	1-245
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.