



Recombinant Mouse Tetraspanin-6 (Tspan6)

Product Code	CSB-CF025164MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O70401
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MASPSRRLQTKPVITCLKSVLLIYTFIFWITGVILLAVGIWGVKVSLENYFSLNEK ATNV PFVLIGTGTVIILLGTFGCFATCRTSAWMLKLYAMFLTLIFLVELVAAIVGFVFRH EIKN SFKSNYENALKEYNSTGDYRSEAVDKIQSTLHCCGVTNYGDWKGVTNYYSETG FPKSCCKL EGCYPQRDADKVN EEGCFIKVMTTIESEMGVVAGISFGVACFQLIGIFLAYCLS RAITNN QYEIV
Source	in vitro E.coli expression system
Target Names	Tspan6
Protein Names	Recommended name: Tetraspanin-6 Short name= Tspan-6 Alternative name(s): Transmembrane 4 superfamily member 6
Expression Region	1-245
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.