



Recombinant Rhizobium meliloti Flagellar biosynthetic protein fliP (fliP)

Product Code	CSB-CF339877RKU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P37827
Product Type	Transmembrane Protein
Immunogen Species	Rhizobium meliloti (strain 1021) (Ensifer meliloti) (Sinorhizobium meliloti)
Sequence	MLRFATFIIAMMAMSGIAGAQSFPADILNTPVDGSVASWIIRTFGLLTVLSVAPGI LIMV TSFPRFVIAFAILRSGMGLATTPSNMIMVSLALFMTFYVMAPTFDRAWRDGIDP LLKNEI SETDAMQRMSEPFREFMVANTRDKDLQLFIDIAREKGQTVVVDEKVDLRAVVP AFMISEI RRGFEIGFLIMLPFLVIDLIVATITMAMGMMMLPPTAISLPFKILFFVLIDGWNLLV GSL VRSFI
Source	in vitro E.coli expression system
Target Names	fliP
Protein Names	Recommended name: Flagellar biosynthetic protein fliP
Expression Region	1-245
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.