



Recombinant Hordeum vulgare ATP synthase subunit a, chloroplastic (atpl)

Product Code	CSB-CF374081HWQ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A1E9I5
Product Type	Transmembrane Protein
Immunogen Species	Hordeum vulgare (Barley)
Sequence	MNIIPCSIKTLKGLYDISGVEVGQHFYWQIGGFQIHAQVLITSWVVITILLGSVVIA VRN PQTIPTDGQNFFFEYVLEFIRDLSKTQIGEEYGPWVPPFIGTMFLFIFVSNWSGALL PWKII ELPHGELAAPTNDINTTVALALLTSAAYFYAGLSKKGLSYFEKYIKPTPILLPINIL EDF TKPLSLSFRLFGNILADELVVVVLVSLVPLVIPIVMFLGLFTSGIQALIFATLAAA YIG ESMEGHH
Source	in vitro E.coli expression system
Target Names	atpl
Protein Names	Recommended name: ATP synthase subunit a, chloroplastic Alternative name(s): ATP synthase F0 sector subunit a F-ATPase subunit IV
Expression Region	1-247
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.