



Recombinant Human Transmembrane protein 41A (TMEM41A)

Product Code	CSB-CF023837HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q96HV5
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	LYLLSTRLPRGRRLGSTEEAGGRSLWFPSDLAELRELVREYRKEHQAYVF LLFCGAY LYKQGFAIPGSSFLNVLGALFGPWLGLLLCCVLTSVGATCCYLLSSIFGKQLV VSYFPD KVALLRKVEENRNSLFFLLFLRLFPMPNWFNLSAPILNIPIVQFFFVLI GLI PYN FICVQTGSILSTLTSLDALFSWDTVFKLLAIAMVALIPGTLIKKFSQKHLQLNETS TANH IHSRKDT
Source	in vitro E.coli expression system
Target Names	TMEM41A
Protein Names	Recommended name: Transmembrane protein 41A
Expression Region	18-264
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.