



Recombinant Phalaenopsis aphrodite subsp. formosana ATP synthase subunit a, chloroplastic (atpl)

Product Code	CSB-CF658956PDAO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q3BAQ4
Product Type	Transmembrane Protein
Immunogen Species	Phalaenopsis aphrodite subsp. formosana (Moth orchid)
Sequence	MNVIPCSIKTLKGLYDISGVEVGQHLYWQIGGLQIHAQVLITSWVVIAILLGSVIIA VRN PQTIPTNGQNFFEYVLEFIQDLSKTQIGEEYGPWVPPFIGTMFLFIFVSNWSGALL PWKII KLPHGELAAPTNDINTTVALALPTSVAYFYAGLRKKGLGYFGKYIQPTPILLPINI LEDF TKPLSLSFRLFGNILADELVVVVLVSLVPSVVPPIVMFLGLFTSGIQALIFATLAAA YIG ESMEGHH
Source	in vitro E.coli expression system
Target Names	atpl
Protein Names	Recommended name: ATP synthase subunit a, chloroplastic Alternative name(s): ATP synthase F0 sector subunit a F-ATPase subunit IV
Expression Region	1-247
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.